

ACHIEVE HIGH SURVIVAL RATES WITH

The Cryotec Method

Cryotec products are designed for safe and effective operability

Plates designed for safer embryos manipulation



Built-in groove for Cryotec to keep it stable


Waste-fluid space around the well for efficient operations

A well designed to prevent bubble formation during the warming process

Round-bottom eliminates shadows caused by light refraction

The round shape of the wells facilitates gentle osmolarity change

Cryotec optimized for ease of operation



A thin sheet for maximum rate of cooling and warming

Easy grip long handles with plenty of space to note down patient's information

The right-angled triangle black tip makes it easier to distinguish between sides and to see the tip of the Cryotec for easy capping

Solutions enable safe recovery

- Free from serum and SSS, and with low impact on cells due to the use of endotoxin-free trehalose and a low concentration of CPA.
- We use Hydroxypropyl cellulose (HPC) to increase the viscosity of our solutions and slow down the movement of oocytes and embryos, reducing the risk of damage and improving the overall safety of the recovery process.



The Cryotec Method is...

a **safe and effective** technique for vitrification and warming of oocytes and embryos. It offers a unified protocol for oocytes and embryos of all stages leading to **high survival rates**.



NO TIME PRESSURE



ONE PROTOCOL FOR ALL



GENTLE ON OOCYTES AND EMBRYOS



IMPROVED OPERABILITY

Our Products

Vitrification Solution Set 210

(For 10 applications)



1.8 ml Equilibration Solution (ES) x 2
1.8 ml Vitrification Solution (VS) x 4

Produced by Reprolife America, LLC

Warming Solutions Set 220

(For 10 applications)



1.8 ml Thawing Solution (TS) x10
1.8 ml Diluent Solution (DS) x 2
1.8 ml Washing Solution (TS) x4

Produced by Reprolife America, LLC

Warming Solutions Set 120

(For 5 applications)



1.8 ml Thawing Solution (TS) x 5
1.8 ml Diluent Solution (DS) x 1
1.8 ml Washing Solution (WS) x 2

Produced by Reprolife America, LLC

Vitri Plate

(Set of 10)



Produced by REPROLIFE Inc.

Warm Plate

(Set of 10)



Produced by REPROLIFE Inc.

Cryotec

10 Cryotec vitrification devices



Produced by REPROLIFE Inc.



Available in blue, green, yellow, white and red



How to order

Please contact Cryotec America, LLC by email orders@cryotecamerica.com

This brochure was updated in August 2025